



denovodental.com

ISO
ISO
Certified

CE
CE
Certified

FDA
FDA
Registered

Autoclavable
Autoclavable

Surgical
Surgical
Stainless Steel

★
Made in
the USA



General Purpose Molar Bands

Application Procedures



1.800.854.7949



1.800.847.8599



info@denovodental.com



Responsibly Printed
& Packaged

Made and assembled in the USA
DENOVO Dental Inc. Baldwin Park, CA 91706




Molar Band Application

1 TRIAL FIT & SEAT

TRIAL FIT: Select an approximate upper or lower band from the provided assortment of sizes and trial fit to tooth. A properly seated band will be difficult to dislodge or push gingivally.

If trial fit is incorrect, **sterilize the band for future use***.



SEAT & CEMENT: The entire band and any attachments should be free of occlusal contact with the opposing arch. Cement band to tooth with your preferred glass ionomer cement.

**Bands can be sterilized using autoclave, chemiclave or dry heat.*

Molar Band Kit Refills


Refilling your kit is easy! Orders can be placed online at www.denovodental.com or over the phone by calling: **(800) 854-7949**.

General Purpose Molar Bands are sold individually to easily refill the sizes that you use most.



Molar Band Refills

TYPE:	SIZE:
MAXILLARY	23 – 42
MANDIBULAR	23 – 42



All Denovo products available to order at www.denovodental.com

Have additional questions? Talk to the experts! Call us at:
1.800.854.7949

FREE Ground Shipping on Orders Over \$300

Continental US only, certain restrictions apply, see website for full details

 denovodental.com



CE Certified
2195



Denovo Dental Inc.
5130 Commerce Dr.
Baldwin Park, CA, USA



KinderDent, Gutenbergstrasse 7
D – 28844 Weyhe, Germany



Denovo Bands are made from surgical grade stainless steel, which contains nickel (9.58%). Nickel has been shown to cause allergic reactions in some patients.



Do not reuse. Appliances are intended for single use only. If appliance is only used for a trial fit, sterilize for future use.



Aspiration/Ingestion possibility. Use rubber dam on installation and educate patient to avoid hard or chewy substances which could dislodge appliance.

Denovo products are sold clean but not sterile; sterilization prior to use is recommended for patients with compromised immune systems.

REV [1]