



## Molar Band Order Form

| Upper (Maxillary) | Size | Lower (Mandibular) |
|-------------------|------|--------------------|
| 146               | 23   | 246                |
| 147               | 23½  | 247                |
| 148               | 24   | 248                |
| 149               | 24½  | 249                |
| 150               | 25   | 250                |
| 151               | 25½  | 251                |
| 152               | 26   | 252                |
| 153               | 26½  | 253                |
| 154               | 27   | 254                |
| 155               | 27½  | 255                |
| 156               | 28   | 256                |
| 157               | 28½  | 257                |
| 158               | 29   | 258                |
| 159               | 29½  | 259                |
| 160               | 30   | 260                |
| 161               | 30½  | 261                |
| 162               | 31   | 262                |
| 163               | 31½  | 263                |
| 164               | 32   | 264                |
| 165               | 32½  | 265                |
| 166               | 33   | 266                |
| 167               | 33½  | 267                |
| 168               | 34   | 268                |
| 169               | 34½  | 269                |
| 170               | 35   | 270                |
| 171               | 35½  | 271                |
| 172               | 36   | 272                |
| 173               | 36½  | 273                |
| 174               | 37   | 274                |
| 175               | 37½  | 275                |
| 176               | 38   | 276                |
| 177               | 38½  | 277                |
| 178               | 39   | 278                |
| 179               | 39½  | 279                |
| 180               | 40   | 280                |
| 181               | 40½  | 281                |
| 182               | 41   | 282                |
| 183               | 41½  | 283                |
| 184               | 42   | 284                |

DENOVO general purpose molar bands are made in regular temper and narrow occlusal- gingival height. These bands, designed with maxillary and mandibular shapes, have universal anatomy to fit both 1st primary, 2nd primary, and 1st permanent molars. Each band is fully proportioned to provide a gradual increase in height as band size increases.

Sizes are available from 23-42, with ½ sizes available for precise fitting. The bands are completely festooned.

To use the order form at left:

1. Enter the quantities for each item in the white boxes (upper bands on the left, lower bands on right).
2. If you require complete part numbers for a purchase order, add the prefix 160 to all part numbers (an upper size 25 band has a part # of 160-150).
3. Quantity discounts are available for your total order. Add up all of the bands ordered to find out

| Total Bands ordered | Price per band |
|---------------------|----------------|
| 1-29 bands          | Retail         |
| 30-49 bands         | 5% discount    |
| 50-99 bands         | 10% discount   |
| 100 or more         | 15% discount   |

### Molar Band Kits

|                                   |         |  |
|-----------------------------------|---------|--|
| Molar Band Kit Upper & Lower      | 161-100 |  |
| Molar Band Kit Maxillary (Upper)  | 161-101 |  |
| Molar Band Kit Mandibular (Lower) | 161-102 |  |

| Kit Description                                      | Part #  |
|--|---------|
| Maxillary & Mandibular Band Kit (Contains 288 Bands) | 161-100 |
| Maxillary Band Kit (Contains 144 Bands)              | 161-101 |
| Mandibular Band Kit (Contains 144 Bands)             | 161-102 |

Name \_\_\_\_\_

Address \_\_\_\_\_

City, State, Zip \_\_\_\_\_

Phone \_\_\_\_\_

Comment \_\_\_\_\_