



### **Composition of DENOVO Stainless Steel 304**

(Used for temporary space maintenance or temporary protection of primary teeth)

Iron	69.22%
Chromium	18.12%
Nickel	9.58%
Manganese	1.59%
Silicon	0.52%
Molybdenum	0.39%
Copper	0.35%
Cobalt	0.18%
Carbon	0.05%

Trace amounts of the following:  
Phosphorus, Titanium, Sulphur, Aluminum, Oxygen.

For purposes of **MRI compatibility**, the material is composed primarily of iron; and are affected by magnetic fields.